

SCS 910- O² CYLINDER SUPPORT



This support was developed with the objective of facilitating the transport of the O² cylinder during the removal of patients who need respiratory assistance.

FEATURES

- Structure in painted carbon steel sheet
- O₂ Velcro belts for fixing the cylinder
- Allows attachment in different positions on the stretcher
- Holds cylinders up to 230mm in diameter
- Includes docking system for ambulance transport
- Compatible with the MRS 310 / MWS 320 / MXS 330 / MCS 360 stretchers

TECHNICAL INFORMATIONS

- Total height: 508mm
- Total width: 139mm
- Total length: 176mm
- Net Weight: 1.8Kg
- Carrying Capacity: 10.00Kg