

MPS 440- HIBERNA®

CAMPAING BED



The MCS 440 Hiberna® camp bed was developed for use in field hospitals, as an additional option for patient accommodation. Hiberna® has a hollow mesh base, which allows for a better fit and fluid flow, in addition to having a folding closing system and a bag for transport and storage.

FEATURES

- Folding structure in hardaluminium
- Base in hollow and removable mesh, which allows the flow of fluids
- Stabilizer bar for opening and closing
- Finishing in e-coat paint

STANDARDS ACCESSORIES

- Carrying bag

TECHNICAL INFORMATIONS

- Total open length: 2000mm
- Total folded length: 1000mm
- Total open width: 735mm
- Total folded width: 230mm
- Open height: 514mm
- Closed height: 100mm
- Net Weight: 7.6Kg
- Carrying Capacity: 220Kg

STANDARDS/ CERTIFICATIONS

- STANAG 2040
- Anvisa Registry nº 80438490002
- AFE Anvisa - 8043849

