

# MTS 410- PEGASUS®

PEGASUS MILITARY TRANSPORTER

**SITMED**  
NUMBER ONE IN LIFE TRANSPORT



The military transport MTS 410 Pegasus®, was developed specifically to attach the Pandora III® multi-folding stretcher, and can be used in TACEVAC (Tactical Evacuation), CASEVAC (Casualty Evacuation) and MEDEVAC (Medical Evacuation) situations, and in the transport of injured in catastrophes of any nature and dimension as well as for employment in field hospitals.

Pegasus® has an opening and closing system, ensuring the versatility and agility necessary for military operations, in addition to having two 20' rim wheels, ideal for use in irregular and uneven terrain.

In addition to carrying the Pandora III® multi-folding stretcher, the Pegasus® carrier also allows the attachment of the Meduse® immobilization board, for situations where a complete immobilization of the patient is necessary.

## FEATURES

- Structure in hard aluminum
- Detachable 20' rim wheels
- Multi-folding system for storage and transport
- Adjustable coupling system, for stretcher or immobilization board
- Approximate assembly time: 10 seconds

## TECHNICAL INFORMATIONS

- Total open length: 1326mm
- Total folded length: 400mm
- Total open width: 696mm
- Total folded width: 680mm
- Total open height: 971mm
- Total folded height: 800mm
- Net Weight: 25.00 Kg
- Carrying Capacity: 200Kg

## STANDARDS/ CERTIFICATIONS

- STANAG 2040
- NATO Stock Number 6530-19-006-9262
- CE certification
- Anvisa Registration nº 80438490002
- AFE Anvisa - 8043849

