


# MPS 120- PANDORA II®

MULTI-FOLDING STRETCHER



The Pandora II® folding stretcher was developed to facilitate the manual transport of victims in the most diverse rescue situations, through a light, practical and compact equipment.

Pandora II® has a foldable structure that allows it to reduce its dimensions for storage, and comes with a carrying bag.

## FEATURES

- Structure in hard tubular aluminum
- Base in resistant, washable and waterproof vinyl canvas
- Multi-folding system for storage
- 04 rubber shoes for support

## STANDARD ACCESSORIES

- 02 Quick-connect PVC safety belts
- Bag with carrying handle

## TECHNICAL INFORMATIONS

- Total open length: 2086mm
- Total folded length: 1053mm
- Total folded width: 130mm
- Total open height: 110mm
- Total folded height: 180mm
- Net Weight: 5.20Kg
- Carrying Capacity: 160.00Kg
- Vinyl Color: Black

## STANDARD / CERTIFICATIONS

- DIN EN 1865
- CE certification
- Anvisa Registration nº 80438490002
- AFE Anvisa - 8043849

