

MRS 310 - ESSENTIAL®

RETRACTABLE AMBULANCE STRECHER



The MRS 310 Essential® stretcher is a retractable operating stretcher that can be custom-made for the vehicle to which it is to be installed. It has a light, reinforced structure, ideal for routine patient rescue and removal. It has a unique set of leg retraction, which retract when activated to enter a vehicle, and automatically arm when taking it out of the ambulance, offering agility, safety and practicality to rescuers.

FEATURES

- Duralumin structure
- Backrest adjustable in 6 positions
- Safety lock that prevents accidental disarming
- 04 swivel castors with brakes
- Swinging side handles
- Retractable folding system

STANDARD ACCESSORIES

- Waterproof mattress D33 Kg/m³
- 01 four-point automotive safety belt for chest and shoulders
- 02 two-point automotive safety belts for abdomen and legs
- 10G certified

TECHNICAL INFORMATION

- Bed length: 1905mm
- Bed width: 505mm
- Total width: 628mm
- Total length: 1980mm ~2100mm ¹
- Net Weight: 33Kg
- Maximum Load Capacity: 300Kg
- Wheel diameter: ø127mm
- Entry height: 530mm ~ 800mm ²
- Side bars length: 612mm
- Side bars height: 215mm
- Mattress color: Black

OPTIONAL ACCESSORIES

- IV Pole (Ref.: SSS 900)
- O² Cylinder Holder (Ref.: SCS 910)
- Support for Monitors and peripherals (Ref.: SMS 940)
- Trendelenburg anti-shock position
- Neonate Incubator transport version available
- Bariatric Patient transport version available

STANDARDS / CERTIFICATIONS

- AMD Standard 004
- DIN EN 1865
- BS EN 1789
- ABNT NBR 14.561
- 10G Crash-Test
- CE Certification (Class I Medical Device according to Dir. 93/42/CEE)
- Registro Anvisa nº 80438490002
- AFE Anvisa - 8043849

¹ The length of the stretcher when in its "closed" position may vary according to the height needed to adapt to the vehicle to which it is intended to be suited.

² The entry height of the stretcher should be considered according to the height of the floor of the ambulance for which it is to be installed.

