


MPS 420- KRISALIS®

TRANSFER BASE



The MPS 420 Krisalis® transfer base was developed for the rescue of wounded people in combat zones and other hostile relief situations.

Krisalis® is made of washable and highly resistant nylon mesh, which allows it to be dragged or carried through its carrying handles, and is also the ideal accessory for transferring patients between different beds.

FEATURES

- 07 carrying handles
- Base in open weave fabric that allows the flow of fluids
- Flexible and easy-to-storage material

TECHNICAL INFORMATIONS

Total length: 2062mm
Total width: 870mm
Net Weight: 1.2Kg
Carrying Capacity: 150.00Kg

STANDARDS / CERTIFICATIONS

STANAG 2040
CE Certification
Anvisa Registration nº 80438490002
AFE Anvisa - 8043849

