


# MPS 400- PANDORA®

MULTI-FOLDING STRETCHER



The MPS 400 Pandora III® multi-folding stretcher was developed to meet the needs of military evacuation and rescue, having an extremely agile opening and closing system, in addition to being compact and light, ideal for combat emergencies. Manufactured in accordance with NATO standards for aeromedical rescue, the Pandora III® stretcher has retractable handles to make it even more compact and can be used in conjunction with the Pegasus® military transporter, which offers greater mobility when transporting the wounded.

## FEATURES

- Structure in aeronautical alloy aluminum
- Base in open weave fabric that allows the flow of fluids
- Multi-folding closing system
- Shoes with sliding support
- Approximate assembly time: 10 seconds

## STANDARD ACCESSORIES

- Backpack for transport
- Set of seat belts

## TECHNICAL INFORMATION

- Total length open stretcher (extended handles): 2270mm
- Total length open stretcher (handles collected): 2010mm
- Total open stretcher width: 584mm
- Open stretcher height: 145mm
- Folded stretcher dimensions: 512mm x 242mm x 155mm
- Net Weight: 8.2Kg
- Load Capacity: 200Kg

## STANDARD/ CERTIFICATIONS

- STANAG 2040
- DIN EN 1865
- NATO Stock Number 6530-19-006-9263
- CE Certification
- Anvisa Registration No. 80438490002
- AFE Anvisa - 8043849

## OPTIONAL ACCESSORIES

- Serum Support
- Backrest

