

# MCS 200 - ÁGILA®

FOLDABLE RESCUE CHAIR



The MCS 200 Ágila® is a chair designed for use in medical rescues and evacuations, being ideal for use in situations where there is not enough space to use a stretcher, such as: elevators, staircases and narrow corridors. It has a light, reinforced structure, ideal for routine patient rescue and removal. It has a folding lock for storage and a locking system for installation in ambulances.

## FEATURES

- Duralumin structure
- Seat and back in polypropylene material
- 04 Castors with brake, 02 swivels and 02 fixed
- Telescopic front handles and folding rear handles
- Foldable closure
- Safety lock that prevents accidental closing
- 01 02-point automotive safety belt

## STANDARD ACCESSORIES

- Waterproof mat D33kg/m³
- 04 retractable automotive seat belts
- 10G certified

## TECHNICAL INFORMATION

- Open length: 715mm
- Closed length: 323mm
- Width: 485mm
- Height: 1013mm<sup>1</sup>
- Net Weight: 11Kg
- Maximum Load Capacity: 160Kg

## OPTIONAL ACCESSORIES

- Armrest
- Velcro belt for the feet

## STANDARDS / CERTIFICATIONS

- DIN EN 1865
- CE Certification (Class I Medical Device according to Dir. 93/42/CEE)
- Anvisa Registration nº 80438490002
- AFE Anvisa - 8043849

