


SMS 940- MED RAC®

SUPPORT FOR MONITOR AND CARDIOVERTER



The Med Rac® monitor and cardioverter support was developed to facilitate the transport of monitors and peripherals during rescue and removal situations, offering greater safety and convenience. It has side straps for fixing the support to the stretcher, in addition to a folding mechanism to facilitate storage when not in use.

FEATURES

- Folding steel structure
- 04 Velcro straps for securing equipment
- 02 side belts for attachment to the stretcher
- Compatible with the MRS 310 / MWS 320 / MXS 330 stretchers

TECHNICAL INFORMATION

- Total length: 597mm
- Total width: 390mm
- Total open height: 379mm
- Total folded height: 70mm
- Net Weight: 7.0Kg
- Carrying Capacity: 30.0 Kg

