

MXS 330 - ELEVEX®

PANTOGRAPHIC STRETCHER FOR AMBULANCE



The MWS 320 Level Up® stretcher is a bi-articulated stretcher, which can be custom-made for the vehicle to which it is intended to be installed. It has a light, reinforced structure, ideal for routine patient rescue and removal. It has an independent set of leg retraction, which retracts when triggered to enter a vehicle, and automatically arm when removing it from the ambulance. The differential of the independent leg collection system is that it allows a significant reduction in the rescuer's effort to carry out the entry and removal of the stretcher from the ambulance, providing greater ergonomics and facilitating the operation, ensuring more safety for the rescuer and patient.

FEATURES

- Steel and duralumin structure
- Backrest adjustable in 6 positions
- Foldable bed for use in confined spaces
- Bed height adjustable in 8 positions
- 04 swivel castors with brakes
- Swing side handles
- Pantographic leg recoil system

STANDARDS ACCESSORIES

- Waterproof mat D33 Kg/m³
- 01 04-point automotive safety belt for chest and shoulders
- 02 02-point automotive safety belts for abdomen and legs
- Bed height adjustable in 8 positions
- 10G certified anchorage system
- Trendelenburg anti-shock positions
- Rigid backing for CPR operations

TECHNICAL INFORMATION

- Bed length: 1905mm
- Total length: 2015mm
- Length with reduced bed: 1570mm
- Bed width: 505mm
- Total width: 635mm
- Net Weight: 56.5Kg
- Maximum Load Capacity: 300Kg
- Wheel diameter: ø200mm
- Entry height: 450mm ~ 880mm
- Length of the side straps: 612mm
- Height of side straps: 215mm
- Mat Color: Black

OPTIONAL ACCESSORIES

- Serum and Blood Support (Ref.: SSS 900)
- Support for O₂ Cylinder (Ref.: SCS 910)
- Support for Monitors (Ref.: SMS 940)
- Trendelenburg anti-shock position
- Additional side anchoring system

STANDARDS / CERTIFICATIONS

- AMD Standard 004
- DIN EN 1865
- BS EN 1789
- ABNT NBR 14.561
- 10G Crash-Test
- CE Certification (Class I Medical Device according to Dir. 93/42/CEE)
- Anvisa Registration nº 80438490002
- AFE Anvisa - 8043849

